GISKIT MD B.V.

Instructions for Use

REF GIS-Kit: GISK-LM969 or GISK-LN970

<p>ExEm® Gel</p> <p>Giskit MD B.V. Veerkade 5 F 3016 DE Rotterdam, The Netherlands</p> <p>CE 2460</p> <p>STERILE ! MD</p>	<p>GIS-Catheter.</p> <p>Gynetics Medical Products N.V. Rembert Dodoensstraat 51 3920 Lommel, Belgium</p> <p>CE 1639</p> <p>STERILE EO MD PHT DEHP</p>
<p>⚠ LATEX ⚠ ⚠ ⚠ ⚠ ⚠ 30°C</p>	
<p>Assembled by Giskit MD B.V., Veerkade 5 F, 3016 DE Rotterdam, The Netherlands tel: +31 10 2619100 • sales@iq-medicalventures.com • www.iq-medicalventures.com</p>	

Report serious incidents immediately to safety@iq-medicalventures.com

Description:

Gel Instillation Sonohysterography (GIS). The GIS-Kit allows for quick and simple filling of the uterine cavity with ExEm® Gel in order to subsequently obtain high quality ultrasound images using a transvaginal ultrasound probe. The GIS Kits are packed as 12 sets in a box with individually sterile packed syringes ExEm® Gel and individually sterile packed GIS-catheters #5165-LM969 or GIS-catheters #5165-LN970.

Indications and usage: All indications for sonohysterography.

The GIS-Kit is intended for use by medical professionals skilled in sonohysterography.

Contraindications:


The GIS-Kit should not be used in the presence of active pelvic infection, sexually transmitted disease, profuse bleeding or pregnancy. Nor in case of allergy to any of the constituents of ExEm® Gel (hydroxyethylcellulose, glycerol and purified water).

Precautions and warnings:

- The GIS-Kit should be considered as contaminated medical waste after use.
- Do not use if package of ExEm® Gel or GIS-catheter has been previously opened or damaged.
- The GIS-Kit is intended for single use only. Do not reuse, re-sterilize or autoclave. Reuse of single use devices creates a potential risk of patient or user infections. Contamination of the device may lead to injury, illness or death of the patient. Cleaning, disinfection and sterilization may compromise essential material and design characteristics leading to device failure. The manufacturer will not be responsible for any direct, incidental or consequential damages resulting from re-sterilization or reuse.
- Do not use the GIS-Kit after the expiry date of the ExEm® Gel, and/or the expiry date of the catheter and/or the expiry date of the cervical applicator.
- The ExEm® Gel is not to be used for injections.

Instructions for use:

1. Connect a syringe filled with ExEm® Gel suitable for sonohysterography, to the purple handle of the inner tubing/catheter (Fig. 1).
2. Fill the cervical applicator with ExEm® Gel until the tubing and applicator are free of air bubbles (Fig. 2).
3. Insert the cervical applicator into the cervical canal of the uterus until a certain resistance is felt and instill the ExEm® Gel into the uterus until the patient experiences a certain feeling of pain (Fig. 3).
4. Disconnect the syringe and leave the cervical applicator in place by pulling the syringe (and connected inner tubing/catheter) and holding the outer tubing of the catheter; some drops of ExEm® Gel may leak (Fig. 4).
5. Complete the ultrasound procedure.
6. Remove the cervical applicator by pulling the cord, this can be done by the ultrasonographer or the patient (Fig. 5).

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